

CRISPR/Cas9 using the Cas9-GFP Cells

1. Cas9-GFP MSCVneo HA-ER-HoxB8 cells
 - a. Derived by transduction of bone marrow from the Cas9-GFP transgenic mouse
 - b. JAX: <https://www.jax.org/strain/026175>
 - c. These cells have constitutive expression of Cas9 and GFP (bicistronic sequence)

Transduction

1. The cells are readily transducible (infectable) using retrovirus or lentivirus.
 - a. I have typically used lentivirus produced by the Broad Institute.
 - b. These are sgRNA in the pXPR_BRD003 backbone, also known as “lentiguide”.
 - c. This is their most common backbone, and it contains a PURO-R cassette.
2. The Broad Institute supplies virus in ~65 ul aliquots.
 - a. Their viral titer generally approximates 10^7 /ml.
 - b. Therefore, I will use between 30-60 ul (~500,000 viral particles) per transduction
3. Transduction can be carried out in many formats.
 - a. Where possible, I like to use 12-well plates.
 - b. However, we also routinely transduce cells in 96-well plate format.

Transduction of Cas9-GFP ER-HoxB8 cells in 12-well plates

1. Count the cells & resuspend at 500,000/ml.
 - a. Add polybrene to a final concentration of 8.8 ug/ml.
 - b. Aliquot 500 ul (250,000 cells) into each well of a 12-well plate.
2. Add 50 ul of virus to each well
 - Final polybrene conc. = 8 ug/ml
3. Spinfection
 - a. 1000g x 60 minutes at RT.
4. Following spinfection, add 3 ml of fresh media to each well.
 - Final volume = ~4 ml and final polybrene concentration = 2 ug/ml
5. At 24-hours, I perform a half-media change to further dilute the polybrene
 - a. Aspirate ~2 ml of media from the top of the wells.
 - b. The cells remain undisturbed on the bottom.
 - c. Add 2 ml of fresh media.
 - d. This dilutes the polybrene out to a non-toxic concentration of ~1 ug/ml.
6. At 36-48 hours, I begin PURO selection.
 - a. The cells are in ~4 ml of media.
 - b. I prepare PURO at 50 ug/ml (10x)
 - c. Add 440 ul of PURO to each well.
 - d. Final PURO concentration ~ 5 ug/ml.
7. The cells should be fully selected within ~48-72 hours.

Transduction of Cas9-GFP ER-HoxB8 cells in 96-well plates

1. Count the cells & resuspend at 1,000,000/ml.
 - a. Add polybrene to a final concentration of 16 ug/ml.
 - b. Aliquot 30 ul (30,000 cells) into each well of a 96-well plate.
2. Add 30 ul of virus to each well
 - Final polybrene conc. = 8 ug/ml
3. Spinfection
 - a. 1000g x 60 minutes at RT.
4. Following spinfection, add 200 of fresh media to each well.
 - Final volume = ~260 ul and final polybrene concentration = 2 ug/ml
5. At 24-hours, I perform a half-media change to further dilute the polybrene
 - a. Using a multi-channel pipet, remove the top 150 ul from each well
 - b. The cells remain undisturbed on the bottom.
 - c. Add 150 ul of fresh media.
 - d. This dilutes the polybrene out to a non-toxic concentration of ~1 ug/ml.
6. At 36-48 hours, I begin PURO selection.
 - a. The cells are in ~250 ul of media.
 - b. I prepare PURO at 50 ug/ml (10x)
 - c. Add 25 ul of PURO to each well.
 - d. Final PURO concentration ~ 5 ug/ml.
7. The cells should be fully selected within ~48-72 hours.

Single cell cloning

1. Single cell cloning can be done by (a) FACS or (b) limiting dilution
 - Single cell cloning allows one to screen a series of clones to identify a knock-out
2. FACS
 - 96 well round bottom plates
 - Add 250 ul of media to each well
 - Using FACS, add 1 cell per well
3. Limiting dilution
 - Count the cells multiple times to get a good concentration of live cells
 - Using serial dilutions, dilute the cells to 12 cells/ml
 - In a 50 ml conical with 50 ml of media, add 600 cells (12 cells/ml)
 - Make a 1:3 dilution (15 ml of cells + 30 ml of media), now 4 cells/ml
 - Make a second 1:3 dilution, now 1.3 cells/ml
 - From each 50 ml conical, plate 250 ul per well using a multi-channel pipet
 - Now you will have three plates (3 cells per well, 1 cell per well, and 0.3 cells per well)
 - This will help account for counting errors, and usually one of the plates will have ~20 clones for expansion
4. The clones will grow out and can be visualized by day 10-12
 - I do not change the media during this 12-day period.
 - Visualize the clones by holding the plate up to the light and seeing the clones on the bottom of each well
 - Usually the cloning efficiency is 10-30%, so you will get ~10-30 clones per plate
 - Using a P200, remove the top half of the media and discard
 - Mix the clone up and down and transfer all of the cells to a 48-well plate with fresh media
 - Move from 48 to 24 to 12 to 6 well plate during the expansion process